

Electrotherapy

Therapic 9200 – TH1312

Therapic 9400 – TH1313

Electrotherapy is mainly used for:

- 1
- Training of the neuromuscular system
- 2
- Modulation of pain
- 3
- Control or reduction of spasticity
- 4
- Transdermal release of medicinal substances
- 5
- Improvement or maintenance of joint mobility
- 6
- Promotion of wound healing
- 7
- Resolution of the edema



Technical Specifications

Monophase	max 70 mA/70 V
Diphase	max 70 mA/70 V
Sync. Monophase	max 70 mA/70 V
Sync. Diphase	max 70 mA/70 V
Short period	max 70 mA/70 V
Long period	max 70 mA/70 V
Rectangular	max 100 mA/100 V
Exponential	max 100 mA/100 V
Triangular	max 100 mA/100 V
Traebert	max 100 mA/100 V
Rectangular faradic	max 80 mA/80 V
Modulated faradic	max 100 mA/100 V
Triangular neodynamic	max 80 mA/80 V
Neodynamic	max 80 mA/80 V
Tens	max 100 mA/100 V
Tens S/A/R	max 100 mA/100 V
Tens random S/A/R	max 100 mA/100 V
Tens burst S/A/R	max 100 mA/100 V
Biphasic S/A	max 100 mA/100 V
Galvanic	max 50 mA/50 V
Iontophoresis	max 50 mA/50 V
Kotz (Russian stimulation)	max 100 mA/100 V
Interferential	max 100 mA/100 V
Classic Interferential	max 100 mA/100 V
Isoplanar Interferential	max 100 mA/100 V
Vectorial Interferential	max 100 mA/100 V



Therapic 9200

Operation	Constant Voltage Constant Current
Peak power	100 V
Peak current	100 mA
Independent output channels	2
Stored protocols	yes
Storable protocols	200
Ultrasound and vacuumtherapy connection	yes
Display	Graphic colour 320 x 240 px Touch & Scroll
Dimensions - Weight	39 x 14 x 30 cm - 4,4 kg
60x85 mm electrodes	4
Sponges for 60x85 mm electrodes	4
50x50 mm electrodes	4
Sponges for 50x50 mm electrodes	4
10 cm elastic band	2
6 cm elastic band	2



Therapic 9400

Operation	Constant Voltage Constant Current
Peak power	100 V
Peak current	100 mA
Independent output channels	4
Stored protocols	yes
Storable protocols	200
Ultrasound and vacuumtherapy connection	yes
Display	Graphic colour 320 x 240 px Touch & Scroll
Dimensions - Weight	39 x 14 x 30 cm - 4,4 kg
60x85 mm electrodes	8
Sponges for 60x85 mm electrodes	8
50x50 mm electrodes	8
Sponges for 50x50 mm electrodes	8
10 cm elastic band	4
6 cm elastic band	4