



RF Micro-Needle

The RF Micro-Needle handpiece applies 2MHz RF energy to the Micro-Needle. Skin improvements such as skin regeneration and lifting by collagen & elastin regeneration.



Vacuum Technology Accurate Treatment


RF energy is precisely transmitted to target depth and treatment area.

Depth Adjustable Smooth Adjustment

Adjustable needle depth of 0.5mm from 0.5-4.5mm.

Various Multi Needles

Insulated, Non-insulated 10pin, 25pin, 35pin
For various skin conditions and treatment areas




HIFU

Skin lifting procedure by immediately contracting the SMAS (facial band) with the formation of a coagulation point of 65-67°C through 4MHz and 7MHz of high intensity focused ultrasound energy. Applied with Transcure Transducer technology for stable energy delivery.



Technology



Various Cartridges

1.5mm, 3.0mm, 4.5mm, 6.0mm, 1.0cm
for various skin layers & treatments



Before & After

