




POWER SHAPE 2

Roll Handpiece

The combination of 2MHz RF, 658nm Low Level Laser & Vacuum technology for body contouring. In particular, the automatic roller and the vacuum function helps to improve cellulite.



Multi-polar Handpiece

Combination of 2MHz Multi-polar RF, Vacuum & LED. In particular, vacuum function allows effective removal of wastes through lymph circulation. Various size of handpieces enable whole body treatment.



Thermal Handpiece


Treatment on narrow and curved areas is possible with 2MHz Bi-polar RF in rubbing method. Two size of tips can be used selectively depending on the preferred treatment area.



Large



Small



Before & After



Age 35, Treatment 10 times