



FULLY AUTOMATIC LIFTER

It has a load capacity of 150 kg and 8 safety levels, making it not only the top of the range, but the universal reference in pool lifts: the best synthesis of quality, innovation, reliability and technological development.



TECHNICAL FEATURES

- 2 speeds motorised movements
- Control bar far forward motion and seat
- Starter key
- Battery charge level indicator
- Emergency stop button
- Anti-crushing system for backward movement
- Automatic braking system integrated in the motorised wheel
- Electrohydraulic box
- Emergency manual pump for use in case of power failure
- Front rubber stoppers far forward safety lock
- Infrared front hole sensor for forward safety stop
- Retractable leg support
- Adjustable armrests
- Internal battery charger with LED indicator
- Average autonomy: 50 working cycles
- 2-point safety belt
- Stainless steel frame for in-water parts
- Galvanization and powder coating

+ TECHNICAL DETAILS

+ CERTIFICATIONS AND PATENT

+ DOWNLOAD THE INFORMATION SHEET IN PDF

Share this product



MOBILE POOL LIFTS



BluOne
Portable pool lift



BluPool
Portable pool lift



PandaPool
Portable pool lift



Pelican Pool
Portable pool lift



EcoPool
Portable pool lift