



RF Micro-Needle

The RF Micro-Needle handpiece applies 2MHz RF energy to the Micro-Needle. Skin improvements such as skin regeneration and lifting by collagen & elastin regeneration.




Vacuum Technology Accurate Treatment

RF energy is precisely transmitted to target depth and treatment area.

Dropping Motor Smooth Adjustment

Adjustable needle depth of 0.1mm from 0.5–3.5mm.



Insulated

Non-Insulated

Various Multi Needles

Insulated/Non-Insulated 10pin, 25pin, 36pin
For various skin conditions and treatment areas



Cryo Handpiece

Pain management using Cryo


Cryo-anesthesia reduce the down time after treatment.



Large Size



Before & After



Age 150s Treatment 12 times