



SAVER ONE SEMI-AUTOMATIC

Categories: Saver One Series, Uncategorized

Description

Technical specifications

Download

User-friendly PAD (public access defibrillator) semi automatic version.

The rescuer must press the shock button to deliver the shock!

SAVER ONE PADs are designed for a public access use and licensed to administer fast and safe rescues. Highly-effective and user-friendly for any lay rescuer, even without minimal training.

Designed to help save lives in any public circumstances and optimized for a rapid delivery of defibrillation therapy, anytime and anywhere and before EMS team arriving.

Practical and very simple to use. Highly intuitive with pictograms and flashing icons and clear voice messages to lead and support ordinary rescuers through each step for an effective lifesaving protocol.

MODEL NUMBERS

Code SVO-B0001:	Semi- Automatic Standard Version at 200J
-----------------	--

Code SVO-B0002:	Semi- Automatic Power Version at 360J
-----------------	---------------------------------------

CONFIGURATION OPTIONS (BOX CONTENT)

Conf-Norm: Standard Basic Configuration	(adult pads, disposable battery, carrying case)
---	---

Conf-Rech: Rechargeable Configuration	(adult pads, accumulator, charger station, carrying case)
---------------------------------------	---